



X26

Electric
Examination Table

SKD-A



Healthward^{uk}



Back-rest

Height Adjustment

TECHNICAL PARAMETERS

External L*W*H	1810*680* (570-860) mm
Back-rest Adjustment	0-75°

TECHNICAL CONFIGURATION

Pillow	1 pc
Motor	1 pc
Mechanical Lever to Adjust Back Rest	1 pc
2" Casters	2 pcs
Hand Controller	1 pc
Paper Roll Holder	1 pc

ELECTRONICS

Voltage	100-240V
Frequency	50/60Hz
Plug	as required (hospital grade)

QUALITY CONTROL

LEATHER

Flame retardant, yellow resistance, mildew resistance, bending resistance, 100,000 times wear resistance

100,000_{times}

SPONGE

Flame retardant, 4 cm thickness, high-density foam, density 30 kg/m³

4_{cm} | 30_{kg/m³}

PAINTING

11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact

0.12_{mm} | 60° | 50_{kg}

METAL MATERIAL

Baosteel from Fortune Global 500
Panasonic Robotic ensure 360° full smooth welding

500 | 360°

ELECTRIC SYSTEM

Taiwan motor, Noise under 50dB, with lifetime 20,000times movement


50_{db} | 20,000_{times}

BACK ADJUST

Back adjust can choose ratchet or motor

RICH PARTS

LEADING BRAND OF
CHINA MEDICAL FURNITURE



1.2-1.5 mm Thickness powder coated steel frame,
150 kg bear loading




2 pcs 2' Hidden castors, 2 pcs 2' universal wheels



Hand controller; keep qualified after 2000 times
use



Engineering ABS material plug, foot pad




Paper roll holder
Foot bar system
Motor to adjust the back section





INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

CHENGDU SAIKANG MEDICAL EQUIPMENT CO., LTD.

24hrs: 0086-18021231901 
www.health-ward.com

SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,
 Jiangsu Province,China.
 TEL: +86-512-58525850
 FAX: +86-512-58520713
 E-mail: inquiry@jssaikang.com

SAIKANG (USA Branch)

Saikang Medical Technology, Inc
 17800 Castleton St Ste 180,
 City of Industry, Ca 91748
 Tel: +1 949648-2776

SAIKANG (UK Branch)

Uk Healthward International Limited
 Owner(s) Address DEPT 111,
 196 High Road Wood Green,
 London, N22 8HH
 United Kingdom
 TEL: +44 2032399738

SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.
 Viraj Complex Unit 3 Mombasa Road,
 P.O Box 14262-00800, Nairobi.
 TEL: +254 733-704 690

SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC
 X22-S03 England Cluster,
 International City, Dubai UAE
 TEL: +20 128 8962497



Company video



General catalogue

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
 As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.