



Healthward uk

Healthier your life

SKH042

FEATURES AND HIGHLIGHTS

Technical parameters

External size (LxWxH)

800×410×(730-1130)mm



Technical configuration

☑ Air Spring for Height Adjustment	1pc
✓ MDF Table Top	1pc

Application

Movable dining table, nursing bedside table, can be used with various hospital beds.



Height Adjustment

The gas spring system help the table cover moving up and down.



Castors

4 castors made by rubber and could revolve around, convenient, wear-resistant, sliding silent.



Table Top

MDF board on top with PVC side covers.



Bottom Frame

The bottom frame is sprayed with carbon steel and the center of gravity is stable.





















Design

Manufacturing

Maintenance



USA

Saikang Medical Technology, Inc. 17800 Castleton St Ste 180, City of Industry, 91748, California. Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd. Dept. 111, 196 High Road Wood Green, N22 8HH, London. Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC Truda 174, 454080 Chelyabinsk. Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd. Viraj Complex Unit 3, Mombasa Road, P.O. Box 14262-00800, Nairobi. Kenya. Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC X22-S03 England Cluster, International City, Dubai. United Arab Emirates. Tel.: +20 128 8962497

Healthward uk



- No.1, 5th Floor, Building 5, No.360, Wuke West 5th Road •
- +86 18021231901 · \$\sum 24/7\$
- export@health-ward.com
- (www.health-ward.com















Company video

General catalogue