



Healthward uk

Healthier your life

Q5k

FEATURES AND HIGHLIGHTS







Technical parameters

External size (LxWxH):	1930x880x600 mm
Mattress platform:	1710x780 mm
Safe working load:	200 kg
Back-rest adjustment:	0-75°(±5°)
Knee-rest adjustment:	0-45° (±5°)

Technical configuration

☑ Linear motor	2pcs
☑ Aluminium alloy foldable side rail	1set
✓ PP bed ends with safe lock	1set
☑ PP bed platform, removeable, washable	5pcs
☑ 5" double side castors	4pcs
☑ IV pole	1set
✓ IV pole prevision	4pcs
✓ Drainage hook	2pcs

QUALITY CONTROL

SAIKANG medical, as a pioneer in the field of medical furniture, always choose branded suppliers strictly to ensure the safety and reliability of raw material. SAIKANG adopt robot welding, automatic epoxy coating. efficient assembly line, dedicated operator, rigorous quality control, only to guarantee SAIKANG's commitment on quality



Design

This Bed uses unique design to meet hospital strict gremfree request. Head & Foot board and platform are all seamless, which supports effective sterilization.



Steels | welding

Baosteel from Fortune Global 500 Panasonic Robotic ensure 360° ful smooth welding.



Plastic | Painting

Environmental plastic can be sterilized in 100°C, tensile upto 30MPa.

11 Process epoxy painting, ASTM testing anti-baterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact.



Electric Motor

Linear motor,with UL /ROHS/EN standard.Noise under 50 dB,with lifetime 20,000 times movement.

Accessories



Platform

Five sections PP Platform with Integrated mattress Retainer, removable for sufficient sterilization.



Castors | Brake padel

 $5^{\rm m}$ Double covered castors with TPR tire, no worn-out after running 30KM, united forming without bolts

Pass Dynamic test: Bearing 120kg run 30KM and pass obstacles 500time.

The double-folded large-area rubber design, horizontal and vertical two-sided rubber, allows the nurse to touch the brakes when the soles of the feet are soft and large-area plastic surface, saving effort and protecting the upper.













Safe working load

200_{kg}



Back-rest adjustment

75°



Knee-rest adjustment

45°















Design

Research

Manufacturing

Maintenance

Service



South America Branch Saikang Medical Techno-

logy, Inc W3400 Corrientes, Argentina TEL: +5493794312231 Av. Senador Vergueiro, 4204

Rudge Ramos – São Bernardo do Campo – SP, Brazil TEL: +55 (11) 2758-0438

East Africa Branch

Saikang Medical Technology Ltd

Mezzanine Floor, Taj Tower, upper hill, nairobi, Kenya TEL: +254714084806

South East Asia Branch

Saikang Medical Technology Ltd

Sala Thammasop Subdistrict Thawi Watthana District Bangkok Thailand TEL: +66 642741046

UK Branch

Uk Healthward International Limited

DEPT 111, 196 High Road Wood Green, London, N22 8HH United Kingdom TEL: +31 (0) 202111106

Russia Branch

LLC 'SAIKANG' 454004, Russian Federation, Chelyabinsk, Universitetskaya Naberezhnaya, Russia TEL: +7 (919)300-00-06

Saikang Medical Supplies (Mid East) LLC Musaffah P.O Box 53200, Abu Dhabi,

Mid East Branch

UAE Tel: +971 585710652

Healthward uk

Chengdu Saikang Medical Equipment Co., Ltd. STOCK CODE: 870098

- No.1, 5th Floor, Building 5, No.360, Wuke West 5th Road
- → +86 18021231901 · ► 24/7
- export@health-ward.com
- www.health-ward.com







Company video

General catalogue