

V8w5c

ELECTRIC BED




Healthward uk

Healthier your life

V8w5c

FEATURES AND HIGHLIGHTS





Technical parameters

| External size (LxW): | 2085x1020 mm |
|------------------------|--------------|
| Mattress platform: | 1925x990 mm |
| Safe working load: | 250 kg |
| Height adjustment: | 455-745 mm |
| Back-rest adjustment: | 0-75° (±5°) |
| Knee-rest adjustment: | 0-46° (±5°) |
| Trendelenburg: | 0-12° |
| Reverse trendelenburg: | 0-12° |

Optional


| ☐ Auto-regression (back section) |
|--|
| ☐ X-ray translucent mattress platform (back section) |
| ☐ Manual CPR |
| ☐ Night light |
| ☐ Emergency stop button |
| ☐ Battery back-up |
| ☐ Weighing system |
| ☐ Nurse controller |
| ☐ Utility shelf |
| ☐ Trapeze bracket/pole |
| |

Technical configuration

| ✓ Linear motor | 4 pcs |
|---|-------|
| ☑ Hand controller | 1 pc |
| ☑ 5'' Covered Castors | 4 pcs |
| ✓ L-shaped foldable side rail | 1 set |
| ☑ PP bed ends with safe lock with Φ65 bumper wheels | 1 set |
| ☑ PP bed platform with integrated mattress retainer | 1 set |
| ☑ IV pole 4 hooks | 1 set |
| ☑ IV pole prevision | 4 pcs |
| ☑ Drainage hook | 2 pcs |

Optional color (MOQ100pcs-color changing)





Accessories



instant and sufficient sterilization.

Platform

Integrated PP bed platform, adopting advanced manufacturing technology. 4-5 section platform, enables independent installation and easy

Enlarge the seat area when back part is raised, providing more comforts for the patients.

Reduce the depression of the sponge at the seat panel.



Castors

5" covered casters. TPR tire no worn out after running 10 km.
Anti-winding hard shell cover, united forming without bolts.
Pass dynamic test: pass obstacles 40000 times, under the standard 200kg load of the whole bed.



LG Chem fresh PP material, Plug

Equipped with the bumpers wheel, diameter 65mm, preventing



Guardrail

L-shaped foldable guardrail (Patent), bended tube, solid installing, 0 shaking, anti-pinch, ensure patient's safety.





Platform integrated mattress retainer at the both sides of the back and leg section, with a boarder of 33mm height.





standard functions.



Safe working load

 250_{kg}



Back-rest adjustment



Knee-rest adjustment



Height adjustment

 $455 - 745_{\mathsf{mm}}$



Reverse trendelenburg



Trendelenburg















Design

Research

Manufacturing

Maintenance

Service



South America Branch Saikang Medical Techno-

logy, Inc W3400 Corrientes, Argentina TEL: +5493794312231 Av. Senador Vergueiro, 4204

– Sala 31 Rudge Ramos – São Bernardo do Campo – SP, Brazil TEL: +55 (11) 2758-0438

East Africa Branch

Saikang Medical Technology Ltd

Mezzanine Floor, Taj Tower, upper hill, nairobi, Kenya TEL: +254714084806

South East Asia Branch

Saikang Medical Technology Ltd

Sala Thammasop Subdistrict Thawi Watthana District Bangkok Thailand TEL: +66 642741046

UK Branch

Uk Healthward International Limited

DEPT 111, 196 High Road Wood Green, London, N22 8HH United Kingdom TEL: +31 (0) 202111106

Russia Branch

LLC 'SAIKANG' 454004, Russian Federation, Chelyabinsk, Universitetskaya Naberezhnaya, Russia TEL: +7 (919)300-00-06

Mid East Branch

Saikang Medical Supplies (Mid East) LLC


Musaffah P.O Box 53200, Abu Dhabi, UAE Tel: +971 585710652

Healthward uk

Chengdu Saikang Medical Equipment Co., Ltd. STOCK CODE: 870098

No.1, 5th Floor, Building 5, No.360, Wuke West 5th Road





www.health-ward.com







Company video

General catalogue