

SKH048 OXYGEN BOTTLE TROLLEY

Registration code: SKM-B





SKH048

FEATURES AND HIGHLIGHTS

Technical parameters

Size (LxWxH)

1435×360×1040 mm

Technical configuration

☑ 6" TPP Rubber Rastor	2pcs
☑ 3'' Covered Castors	2pcs

Design

U-shaped groove design, the radian meets the requirements of the oxygen bottle for safe transportation.





Size: 1435×360×1040mm, Can load 40L oxygen cylinders, the body is light and easy to move.

Handle

Pushing handle, easy to move the trolley.

6" TPP Rubber Rastor

6" TPP rubber castor, Strong, durable, labor-saving, safty in use.







Steels | Welding Baosteel from Fortune Global 500. Panasonic Robotic ensure 360° full smooth wedlding.



Galvanized Chain 4mm galvanized chain, operational safety.



Thickening Plate Cold rolled steel plate, effective compression, bearing weight, not easy to damage, long service life.

3" Covered Castors 3" Covered Castors, each with brakes, easy to install, stable and smooth.













Design



F



Manufacturing



Maintenance



Service



USA

Saikang Medical Technology, Inc. 17800 Castleton St Ste 180, City of Industry, 91748, California. Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd. Dept. 111, 196 High Road Wood Green, N22 8HH, London. Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC Truda 174, 454080 Chelyabinsk. Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd. Viraj Complex Unit 3, Mombasa Road, P.O. Box 14262-00800, Nairobi. Kenya. Tel.: +254733 704 690

Middle East

Saikang Medical Supplies LLC X22-S03 England Cluster, International City, Dubai. United Arab Emirates. Tel.: +20 128 8962497

Healthward^{`uk}

Chengdu Saikang Medical Equipment Co., Ltd. STOCK CODE: 870098

No.1, 5th Floor, Building 5, No.360, Wuke West 5th Road · Wuhou District (Chengdu) · China

🔲 +86 18021231901 · 🕓 24/7

- ∝ export@health-ward.com
- () www.health-ward.com



