




# SKB041

Operation  
Connecting  
Trolley



Healthward<sup>uk</sup>



Back-rest

Height Adjustment

## MAIN TECHNICAL PARAMETERS

Size	3650*640*(640-970)mm
Back-rest	0-75°(±10°)

## TECHNICAL CONFIGURATION

Crank System	2pc
Back Rest Air Spring	1pc
ABS Side Rail	1set
Directional Wheel	2pc
6" Castors with Central Locking	4pcs
IV Previsions	2pc
IV Pole	1set
Central Locking Break Pedal	4pcs
Directional Wheel Release Pedal	4pcs


### PROTECT

cover all metal screws with rubber cap, protect screws from possible corrosion



### CENTER WHEEL

Can be engaged to facilitate direction control, enabling users to navigate the stretcher easily in hectic and dynamic environments.



### CONNECT WITH OPERATING ROOM

Transfer patient to operation room, prevent cross infection



## QUALITY CONTROL

### METAL MATERIAL

Baosteel from Fortune Global 500  
Panasonic Robotic ensure 360°  
smooth wedding



### PLASTIC

Environmental plastic can be sterilized in 100°C,  
tensile upto 30MPa




### STEELS | PAINTING

1 | Process epoxy painting,  
ASTM testing anti-bacterial, paint  
thickness 0.12mm, brightness  
60°, paint can resist 50kg impact



### METAL DUST COVER

We use metal dust protection for crank system, more durable than plastic cover



### BEARING

We use robust tube in weight bearing structure



### AUTOLOCK

The patient stretcher will be autolocked when it is fixed to one side of the trolley



## RICH PARTS

LEADING BRAND OF  
CHINA MEDICAL FURNITURE



Adjustable vertical oxygen bottle, load capacity: 15KG



Stainless steel handle, more durable



Four double side casters of 6 inch diameter; simultaneously lockable by the step pedal at both sides.


TPR tire no worn out after running 30KM, save anti-winding hard shell, united forming without bolts



Full SS IV pole, bearing 15kg



Water proof mattress with 2 belts





## INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

CHENGDU SAIKANG MEDICAL EQUIPMENT CO., LTD.

24hrs: 0086-18021231901   
[www.health-ward.com](http://www.health-ward.com)

### SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,  
 Jiangsu Province,China.  
 TEL: +86-512-58525850  
 FAX: +86-512-58520713  
 E-mail: [inquiry@jssaikang.com](mailto:inquiry@jssaikang.com)

### SAIKANG (USA Branch)

Saikang Medical Technology, Inc  
 17800 Castleton St Ste 180,  
 City of Industry, Ca 91748  
 Tel: +1 949648-2776

### SAIKANG (UK Branch)

Uk Healthward International Limited  
 Owner(s) Address DEPT 111,  
 196 High Road Wood Green,  
 London, N22 8HH  
 United Kingdom  
 TEL: +44 2032399738

### SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.  
 Viraj Complex Unit 3 Mombasa Road,  
 P.O Box 14262-00800, Nairobi.  
 TEL: +254 733-704 690

### SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC  
 X22-S03 England Cluster,  
 International City, Dubai UAE  
 TEL: +20 128 8962497



Company video



General catalogue

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.  
 As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.