ISO :: D ( E FDA 1/2

# SK-PQ002 MEDICAL BRIDGE CEILING PENDANT (INTEGRATION OF DRY AND WET AREAS)



# Healthward<sup><sup>1</sup>uk</sup>

# Care you precisely

# SK-PQ002

# FEATURES AND HIGHLIGHTS

The medical bridge ceiling pendant (Integration of dry and wet areas) is based on the separation of dry and wet, the two areas are combined into one, suitable for NICU and small nursing areas. Product Application: Suitable for ICU, NICU, PICU, EICU room.

## **Optional Accessories:**

- \* Net basket Stainless steel net basket
- \* We manufacturer and provide different standard of medical gas outlet: British Standard, American standard, German standard, French standard, Japanese Standard, etc.
- \* Different standard of medical gas outlets & adapters: British standard, French standard, German standard, Diss, Ohmeda, Chemetron, and etc.
- \* Inhalers All kinds of oxygen inhalers (disposable, simple, external hanging, etc.)
- \* Negative pressure regulating valve
- \* Waste liquid collection device (external hanging type)
- \* The medical pendant is equipped with individual IV pole, bulit in interface of strong power& weak electronic.
- (Optional cantilever working lamp, monitor holder and surveillance camera in the place of IV pole)

## **Material and Process**

#### **Appearance Design**

The design of bridge hanging style saves ground space; the whole bridge pendant is fully sealed, no sharp angle to prevent the impact.



Imported plastic powder,

electrostatic painting process,

brightness 60°, paint can resist 50kg

painting thickness 0.12mm,

#### Material

Painting

impact.

The tensile strength of aluminum alloy frame is more than 295MPa, some connecting parts are made of cast iron ≥200MPa, and the net load is more than 300kg; one-off extrusion molding, and the surface undergoes the primary oxidation treatment.



Welding Panasonic Robot ensure 360° full smooth welding.



### Function

#### Integration of Dry and Wet Areas

One medical pendant is equipped with medical gas. One part of the medical gas is located in the upper and lower areas, the other part of the medical gas is located in the left and right areas. The design is suitable for small nursing areas with two beds.



#### Electric Translation System

The function column built-in Jiecang motor of model JC35D1, with SAA/TUV/PSE/CE/ UL/KC certificates, ensure safety. Motor noise under 48dB, lifetime 20,000 times, electric translation distance 300mm; the horizontal translation in dry and wet areas, accurate positioning without drift, meeting requirements of different nursing.



#### **Suspension System**

Medical pendant rotation angle 340°, built-in rotation limit device to prevent drift.



## Configuration

#### **Bridge Configuration**

Two fixed LED lightings, 1 set of double cantilever folding working lighting; one power supply and one gas total valve.

Utility Model Patent Negative pressure terminal obtained the Utility Model Patent CN202478404U.

#### 专利内容由知识产权出版社提供



請實:本实用新型涉及一种负压吸引門源终端,包括兩体產盛,在兩体產盛上還定有定位損益。 5所送应位增金上還定有色活圈,在所含色活圈,在所含色活用和國家,在所含品和國家,在所含色活圈,在所含色活圈,在所含色活圈,在所含色活圈,在所含色活用和國家,在所含色活用和國家,在所含色活用和國家,在所含色活用和國家,在所含色活用和國家,在所含色活用和國家,在所含色活用和國家,在所有自己的國家,自己的意思,自己

#### Drawer

Silent self-priming track, no damage after 10,000 times push-pull tests, and the handle position is integrally formed by pressing with mould.



#### **Instrument Platform**

IV Pole with Four Hooks

single hook can bear 4 kg.

The length of IV pole is 600mm, can

be retractable, with four hooks,

The size of the table top is 450\*500mm, and the height of table surface can be adjusted. The table top is made of aluminum alloy panel with iron square tube inside, can bear≥50KG; with external aluminum alloy hanging frame.



**Equipotential Socket** 

Reduce the potential difference between indirect contact voltage and different metal parts.



#### **Medical Gas**

It is equipped with medical gas on both sides of the cabinet and column. The following configurations can be customized, 1. The number of medical gas 2. Interface standard and position Gas port adopts secondary seal, three states (on, off, pull), support maintenance operations while it is working.









Design



F



Manufacturing



Maintenance



Service



USA

Saikang Medical Technology, Inc. 17800 Castleton St Ste 180, City of Industry, 91748, California. Tel.: +1 949648 2776

### United Kingdom

UK Healthward International Ltd. Dept. 111, 196 High Road Wood Green, N22 8HH, London. Tel.: +44 2032399738

#### Russia

Saikang Medical Russia LLC Truda 174, 454080 Chelyabinsk. Tel.: +7 (351) 214 5566 / 8811 / 1516

#### East Africa

Saikang Medical Solutions Ltd. Viraj Complex Unit 3, Mombasa Road, P.O. Box 14262-00800, Nairobi. Kenya. Tel: +2254733 704 690

### Middle East

Saikang Medical Supplies LLC X22-S03 England Cluster, International City, Dubai. United Arab Emirates. Tel.: +20 128 8962497

# **Healthward**<sup>`uk</sup>

Chengdu Saikang Medical Equipment Co., Ltd. STOCK CODE: 870098

#### No.1, 5th Floor, Building 5, No.360, Wuke West 5th Road · Wuhou District (Chengdu) · China

🔲 +86 18021231901 · 🕓 24/7

- ∝ export@health-ward.com
- () www.health-ward.com







General catalogue