

K2k


TWO FUNCTIONS MANUAL BED

Registration code: SK-C1



K2k

FEATURES AND HIGHLIGHTS



Technical parameters

▪ External size (LxW):	2115x980 mm
▪ Mattress platform:	1925x820 mm
▪ Safe working load:	200 kg
▪ Height adjustment:	490mm
▪ Back-rest adjustment:	0-75° (±10°)
▪ Knee-rest adjustment:	0-40° (±10°)

Technical configuration

<input checked="" type="checkbox"/> Manual Cranks	2pcs
<input checked="" type="checkbox"/> 5" Covered Castors	4pcs
<input checked="" type="checkbox"/> Foldable Side Rail	1set
<input checked="" type="checkbox"/> ABS Bed Ends with Safe Lock	1set
<input checked="" type="checkbox"/> IV Pole	1set
<input checked="" type="checkbox"/> IV Pole Prevision	4pcs
<input checked="" type="checkbox"/> Drainage Hook	4pcs
<input checked="" type="checkbox"/> Mattress Holder	2pcs

Optional color (MOQ100pcs- color changing)



QUALITY CONTROL



Design

The design is both modern and practical. The head and foot boards feature a beautiful and elegant shape. The use of a hollow protective wrap angle, combined with a barrel-type bumper, makes the product instantly recognizable.



Steels | welding

Baosteel from Fortune Global 500. Panasonic Robotic ensure 360° full smooth welding.



Plastic | Painting

Environmental plastic can be sterilized at a temperature of 100°C and has a tensile strength of up to 30MPa. 11 Process epoxy painting with a thickness of 0.12mm and a brightness of 60°, which is capable of resisting impacts of up to 50kg.

Accessories



Head & foot board

The head and foot boards feature "integrated protective corners" with hollow design, which are combined with barrel-type bumper wheels that have a diameter of 60mm.

This design not only minimizes the overall space occupied by the bed but also provides protection.

Additionally, the head and foot boards are effortless to disassemble and clean.



Castors

5" covered casters.

PR material, walks 30KM with no wear and tear. Precision-designed anti-winding shell with no screws.

Pass Dynamic test: withstand a weight of 120kg, walk 30KM, and pass through 500 obstacles without any deformation or damage.

Platform

Die punching bed board features a robot-welded surface that has been polished to perfection, leaving no welding scars. Bearing 200kg.

Siderail

The guardrail is made of six-pole folding aluminum alloy, with stainless steel poles and a plastic handle.

Hand Cranks

Two retractable hand cranks to control backrest and kneerest. The hand cranks remain intact after 100,000 uses.

Mattress Retainer

The dual-function mattress retainer serves to not only prevent scratches on the bed frame's paint but also retains the mattress to its place, and prevents mattress sliding.

Drainage Hooks

Injection molded drainage hooks on both sides of the bed.

IV Pole

The infusion pole is made of stainless steel and has four plastic hooks that can bear up to 15 kg of weight.



Safe working load

200 kg



Back-rest adjustment

75°



Knee-rest adjustment

40°



CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

Healthward^{uk}

Chengdu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No.1, 5th Floor, Building 5, No.360, Wuke West 5th Road ·
Wuhou District (Chengdu) · China

☎ +86 18021231901 · 🌐 24/7

✉ export@health-ward.com

🌐 www.health-ward.com



Company video



General catalogue