



Healthward uk

Healthier your life

B6e

FEATURES AND HIGHLIGHTS











Technical parameters

External size (LxW):	2170x1060 mm
Mattress platform:	1925x900 mm
Safe working load:	250 kg
 Height adjustment: 	435-745 mm
■ Back-rest adjustment:	0-75° (±10°)
Knee-rest adjustment:	0-35° (±10°)
Angle indicator:	0-90°
Voltage and frequency:	100-240 V · 50/60 Hz
 Plug as required (hospital grade) 	



Technical configuration

✓ Linear Motor	3 pcs
☑ Hand Controller	1 pc
☑ 5" Hardshell Castors	4 pcs
☑ ABS Side Rail Extra Height	1 set
☑ ABS Bed Ends with Safe Lock	1 set
☑ IV Pole	1 set
✓ IV Pole Prevision	4 pcs
☑ Drainage Hook	2 pcs




Optional

☐ Manual CPR
☐ Electric CPR
☐ Back-up Battery
☐ X-ray Bed Surface
☐ Monkey Pole
☐ Balkan Traction Frame
☐ Directional Wheel



Design

E-shaped large guardrail, with angled display and built-in controller, meets safety standards, providing conforting feeling for patients.



Steels | Painting

Baosteel from Fortune Global 500. 11 Process epoxy painting, ASTM testing anti-baterial, paint thickness 0.12 mm, brightness 60°, paint can resist 50 kg impact.



Welding | Chassis

Panasonic robotic ensure, 360° full smooth wedlding. Solid bed base made of 5.0 mm stamped rotating plate.



Electric Motor

Linear motor, with UL/ROHS/EN standard. Noise under 50 dB, with lifetime 20,000 times movement.



Virgin ABS

Environmental plastic can be sterilized in 100°C, tensile upto 30MPa.



Angle Indicator

Patented angle indicator.



Bed Surface

Stamping forming platform, smooth surface no welding scar, bearing 250 kg load. Option ABS platform and resin X-ray translucent platform.



Side Rail

Full ABS material, tuck away side rail with built-in controller.
From top of the guardrail to bed surface is about 345mm, more height provides better protection for patient. Buttom can endure 100,000times press, buffer, built-in angle indicator.



Head & Foot Board

Full ABS material, Buttom can endure 100,000 times press Head and foot board can be reversable.



Wheels

TPR tire no worn out after running 30KM, save anti-winding hard shell, united forming without bolts.
Pass Dynamic test: Bearing 120kg run 30KM and pass obstacles 500time.



IV pole

Stainless Steel Column with plastic hook.



Drainage Hooks

Mobile plastic drainage hooks.

















Design

Research

Manufacturing

Maintenance

Service



USA

Saikang Medical Technology, Inc. 17800 Castleton St Ste 180, City of Industry, 91748, California. Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC Truda 174, 454080 Chelyabinsk. Tel.: +7 (351) 2145566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd. Viraj Complex Unit 3, Mombasa Road, P.O. Box 14262-00800, Nairobi. Kenya. Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC X22-S03 England Cluster, International City, Dubai. United Arab Emirates. Tel.: +20 128 8962497

Healthward uk



- No.1, 5th Floor, Building 5, No.360, Wuke West 5th Road
- +86 18021231901 · \$\sum_2\$ 24/7
- export@health-ward.com
- www.health-ward.com







Company video

General catalogue