



Healthward uk

Healthier your life

V6k

FEATURES AND HIGHLIGHTS









Technical parameters

External size (LxW):	2085x995 mm
Mattress platform:	1925x900 mm
Safe working load:	200 kg
Height adjustment:	420-730mm
Back-rest adjustment:	0-75°(±5°)
Knee-rest adjustment:	0-46° (±5°)

Optional

☐ Battery Back-up
☐ Bed Extension
☐ Nurse Controller
☐ Emergency Stop Button
☐ Auto-regression Function (Back Section)
☐ Night Light
☐ Manual CPR
☐ X-ray Translucent Mattress Platform (Back Section)
☐ Utility Shelf
☐ Trapeze Bracket/Pole

Technical configuration

✓ Linear motor	3 pcs
☑ Hand controller	1 pc
☑ 5" covered castors	4 pcs
☑ Foldable side rail	1 set
☑ PP bed ends with safe lock with Φ65 bumper wheels	1 set
☑ PP bed platform with integrated mattress retainer	1 set
☑ IV pole 4 hooks	1 set
✓ IV pole prevision	4 pcs
☑ Drainage hook	2 pcs

Optional color (MOQ100pcs-color changing)



QUALITY CONTROL

SAIKANG medical, as a pioneer in the field of medical furniture, always choose branded suppliers strictly to ensure the safety and reliability of raw material. SAIKANG adopt robot welding, automatic epoxy coating. efficient assembly line, dedicated operator, rigorous quality control, only to guarantee SAIKANG's commitment on quality



Design

The unique germ-free design meets hospitals' strict infection control requirment. Seamless PP head and foot board and platform ensure the instant and sufficient sterilization.



Steels | welding

Baosteel from Fortune Global 500. Panasonic Robotic ensure 360° ful smooth welding.



Plastic | Painting

Environmental plastic can be sterilized in 100°C, tensile upto 30MPa.

11 Process epoxy painting, ASTM testing anti-baterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact.



Electric Motor

American motor, with UL/ROHS/EN standard. Noise under 50dB, with lifetime 20,000 times movement.

Accessories



Platform

Integrated PP bed platform, adopting advanced manufacturing technology. 4-5 section platform, enables independent installation and easy maintenance.

Enlarge the seat area when back part is raised, providing more comforts for the patients.

Reduce the depression of the sponge at the seat panel.



Castors

 $5^{\prime\prime}$ covered casters. TPR tire no worn out after running 10 km.

Anti-winding hard shell cover, united forming without bolts.

Pass dynamic test: pass obstacles 40000 times, under the standard 200kg load of the whole bed.















Back-rest adjustment

75°



Knee-rest adjustment

46°



Height adjustment

420-730_{mm}















Design

Research

Manufacturing

Maintenance

Service



South America Branch

Saikang Medical Technology, Inc

W3400 Corrientes, Argentina TEL: +55 (11) 2758-0438 Rua Universal, Hollywood, São Bernardo do Campo, Brazil

TEL: +55 (11) 2758-0438

East Africa Branch

Saikang Medical Technology Ltd

Mezzanine Floor, Taj Tower, upper hill, nairobi, Kenya TEL: +254714084806

South East Asia Branch

Saikang Medical

Technology Ltd
Sala Thammasop Subdistrict
Thawi
Watthana District
Bangkok Thailand
TEL: +66 642741046

UK Branch

Uk Healthward International Limited

DEPT 111, 196 High Road Wood Green, London, N22 8HH United Kingdom TEL: +31 (0) 202111106

Russia Branch

LLC'SAIKANG'

454004, Russian Federation, Chelyabinsk, Universitetskaya Naberezhnaya, Russia TEL: +7 (919)300-00-06

Mid East Branch

Saikang Medical Supplies

(Mid East) LLC Musaffah P.O Box 53200, Abu Dhabi, UAE Tel: +971 585710652

Healthward uk

Chengdu Saikang Medical Equipment Co., Ltd







www.health-ward.com







Company video

General catalogue