



V2k

MANUAL BED

Registration code: SK-C1




Healthward^{uk}

Healthier
your life

V2k

FEATURES AND HIGHLIGHTS



Technical parameters

▪ External size (LxW):	2130x995 mm
▪ Mattress platform:	1925x900 mm
▪ Safe working load:	200 kg
▪ Height adjustment:	500mm
▪ Back-rest adjustment:	0-70° (±5°)
▪ Knee-rest adjustment:	0-38° (±5°)

Technical configuration

<input checked="" type="checkbox"/> Manual crank	2 pcs
<input checked="" type="checkbox"/> 5" covered castors	4 pcs
<input checked="" type="checkbox"/> Foldable side rail	1 set
<input checked="" type="checkbox"/> PP bed ends with safe lock with Ø65 bumper wheels	1 set
<input checked="" type="checkbox"/> PP bed platform with integrated mattress retainer	1 set
<input checked="" type="checkbox"/> IV pole 4 hooks	1 set
<input checked="" type="checkbox"/> IV pole prevision	4 pcs
<input checked="" type="checkbox"/> Drainage hook	2 pcs

Optional

- Auto-regression function (back section)
- Utility shelf
- Trapeze bracket/pole

Optional color (MOQ100pcs- color changing)



QUALITY CONTROL

SAIKANG medical, as a pioneer in the field of medical furniture, always choose branded suppliers strictly to ensure the safety and reliability of raw material. SAIKANG adopt robot welding, automatic epoxy coating, efficient assembly line, dedicated operator, rigorous quality control, only to guarantee SAIKANG's commitment on quality



Design

The unique germ-free design meets hospitals' strict infection control requirement. Seamless PP head and foot board and platform ensure the instant and sufficient sterilization.



Steels | welding

Baosteel from Fortune Global 500. Panasonic Robotic ensure 360° full smooth welding.



Plastic | Painting

Environmental plastic can be sterilized in 100°C, tensile upto 30MPa.
11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact.

Accessories



Platform

Integrated PP bed platform, adopting advanced manufacturing technology.

4-5 section platform, enables independent installation and easy maintenance.

Enlarge the seat area when back part is raised, providing more comforts for the patients.

Reduce the depression of the sponge at the seat panel.



Castors

5" covered casters. TPR tire no worn out after running 10 km.

Anti-winding hard shell cover, united forming without bolts.

Pass dynamic test: pass obstacles 40000 times, under the standard 200kg load of the whole bed.

Head & Foot Board


LG Chem fresh PP material, Plug in type.
Equipped with the bumpers wheel, diameter 65mm, preventing collision.

**Mattress Retainer**


Platform integrated mattress retainer at the both sides of the back and leg section, with a boarder of 33mm height.

**Guardrail**

Aluminum alloy side rail, 6-staff, collapsed type with safe lock.

**Wall bumper wheel**


Horizontal corner bumpers at the outermost periphery of the bed to protect the bed and wall.

**Drainage Hooks**

Mobile plastic drainage hooks.

**IV Pole**

SS IV pole, bearing 15kg.



Safe working load

200 kg



Back-rest adjustment

70°



Knee-rest adjustment

38°



CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service

WORLDWIDE BRANCHES

South America Branch

Saikang Medical Technology, Inc

W3400 Corrientes, Argentina
TEL: +55 (11) 2758-0438
Rua Universal, Hollywood,
São Bernardo do Campo,
Brazil
TEL: +55 (11) 2758-0438

East Africa Branch

Saikang Medical Technology Ltd

Mezzanine Floor, Taj Tower,
upper hill, nairobi, Kenya
TEL: +254714084806

South East Asia Branch

Saikang Medical Technology Ltd

Sala Thammasop Subdistrict
Thawi
Watthana District
Bangkok Thailand
TEL: +66 642741046

UK Branch

Uk Healthward International Limited

DEPT 111,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +31 (0) 202111106

Russia Branch

LLC 'SAIKANG'

454004, Russian Federation,
Chelyabinsk,
Universitetskaya
Naberezhnaya, Russia
TEL: +7 (919)300-00-06


Mid East Branch


Saikang Medical Supplies (Mid East) LLC

Musaffah
P.O Box 53200, Abu Dhabi,
UAE
Tel: +971 585710652

Healthward^{uk}

Chengdu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

 No.1, 5th Floor, Building 5, No.360, Wuke West 5th Road · Wuhou District (Chengdu) · China

 +86 18021231901 ·  24/7

 export@health-ward.com

 www.health-ward.com



Company video



General catalogue