



# V3v5c

## MANUAL BED

Registration code: SK-CD1





**Healthward**<sub>uk</sub>

Healthier  
your life

# V3v5c

## FEATURES AND HIGHLIGHTS



### Technical parameters

External size (LxW):	2130x1025 mm
Mattress platform:	1925x900 mm
Safe working load:	250 kg
Height adjustment:	430-740 mm
Back-rest adjustment:	0-73° (±5°)
Knee-rest adjustment:	0-38° (±5°)
Angle indicator:	0-90°
Auto-regression:	100 mm
Position line:	30°

### Technical configuration

<input checked="" type="checkbox"/> Manual Crank	3 pcs
<input checked="" type="checkbox"/> 5" Covered Castors	4 pcs
<input checked="" type="checkbox"/> PP side rail extra height	1 set
<input checked="" type="checkbox"/> PP bed ends with safe lock with Ø65 bumper wheels	1 set
<input checked="" type="checkbox"/> PP bed platform with integrated mattress retainer	1 set
<input checked="" type="checkbox"/> IV pole 4 hooks	1 set
<input checked="" type="checkbox"/> IV pole prevision	4 pcs
<input checked="" type="checkbox"/> Drainage hook	2 pcs

### Optional

- Utility shelf
- Trapeze bracket/pole

### Optional color (MOQ100pcs- color changing)



# QUALITY CONTROL

SAIKANG medical, as a pioneer in the field of medical furniture, always choose branded suppliers strictly to ensure the safety and reliability of raw material. SAIKANG adopt robot welding, automatic epoxy coating, efficient assembly line, dedicated operator, rigorous quality control, only to guarantee SAIKANG's commitment on quality



## Design

The unique germ-free design meets hospitals' strict infection control requirement. Seamless PP head and foot board and platform ensure the instant and sufficient sterilization.




## Steels | welding

Baosteel from Fortune Global 500. Panasonic Robotic ensure 360° full smooth welding.



## Plastic | Painting

Environmental plastic can be sterilized in 100°C, tensile upto 30MPa.  
11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12mm, brightness 60°, paint can resist 50kg impact.



## Auto-Regression

Auto-Regression 10cm. When repositioning the backrest and thigh rest, the V bed makes the back plate retreat 10cm, so to minimises the pressure and shear forces in the patient's pelvic region. This results in a significant reduction in the development of pressure ulcers.

## Accessories



### Platform

Integrated PP bed platform, adopting advanced manufacturing technology. 4-5 section platform, enables independent installation and easy maintenance.

Enlarge the seat area when back part is raised, providing more comforts for the patients.

Reduce the depression of the sponge at the seat panel.



### Castors

5" covered casters. TPR tire no worn out after running 10 km.

Anti-winding hard shell cover, united forming without bolts.

Pass dynamic test: pass obstacles 40000 times, under the standard 200kg load of the whole bed.

**Head & Foot Board**

LG Chem fresh PP material, Plug in type. Equipped with the bumpers wheel, diameter 65mm, preventing collision.




**Angle Indicator**

With built-in angle indicator, it clearly show back-rest and trendelenburg position.




**Guardrail**

Tuck away side rail with buffer integrated. Integrated handle facilitates the patients to get up/off the bed.



**Wall bumper wheel**

Horizontal corner bumpers at the outermost periphery of the bed to protect the bed and wall.




**Mattress Retainer**

Platform integrated mattress retainer at the both sides of the back and leg section, with a boarder of 33mm height.



**Drainage Hooks**

Plastic drainage hook integrated inside the siderails.




**Drainage Hooks**

Mobile plastic drainage hooks.



**IV Pole**

SS IV pole, bearing 15kg.



Back-rest adjustment

73°



Knee-rest adjustment

38°



Height adjustment

430-740mm



## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service

## WORLDWIDE BRANCHES

### South America Branch

**Saikang Medical Technology, Inc**  
 W3400 Corrientes, Argentina  
 TEL: +5493794312231  
 Av. Senador Vergueiro, 4204  
 - sala 31  
 Rudge Ramos - São Bernardo do Campo - SP, Brazil  
 TEL: +55 (11) 2758-0438

### East Africa Branch

**Saikang Medical Technology Ltd**  
 Mezzanine Floor, Taj Tower,  
 upper hill, nairobi, Kenya  
 TEL: +254714084806

### South East Asia Branch

**Saikang Medical Technology Ltd**  
 Sala Thammassop Subdistrict  
 Thawi  
 Watthana District  
 Bangkok Thailand  
 TEL: +66 642741046

### UK Branch

**Uk Healthward International Limited**  
 DEPT 111,  
 196 High Road Wood Green,  
 London, N22 8HH  
 United Kingdom  
 TEL: +31 (0) 202111106

### Russia Branch


**LLC 'SAIKANG'**  
 454004, Russian Federation,  
 Chelyabinsk,  
 Universitetskaya Naberezhnaya,  
 Russia  
 TEL: +7 (919)300-00-06


### Mid East Branch

**Saikang Medical Supplies (Mid East) LLC**  
 Musaffah  
 P.O Box 53200, Abu Dhabi,  
 UAE  
 Tel: +971 585710652

# Healthward<sup>uk</sup>

Chengdu Saikang Medical Equipment Co., Ltd.  
 STOCK CODE: 870098

 No.1, 5th Floor, Building 5, No.360, Wuke West 5th Road · Wuhou District (Chengdu) · China

 +86 18021231901 ·  24/7

 export@health-ward.com

 www.health-ward.com



Company video



General catalogue