



SKB041-2

PATIENT TRANSPORTATION TROLLEY

Registration code: SKM-F



Healthward^{uk}

Healthier
your life

SKB041-2


FEATURES AND HIGHLIGHTS

Technical parameters

▪ External size (LxW):	1930x640 mm
▪ Table height adjustment:	630-930 mm
▪ Back plate adjustment:	0-75° (±10°)
▪ Reverse Trendelenburg:	0-15°
▪ Trendelenburg:	0-15°

Technical configuration

<input checked="" type="checkbox"/> Hydraulic pump	2 pcs
<input checked="" type="checkbox"/> Back rest air spring	1 pc
<input checked="" type="checkbox"/> ABS side rail	1 set
<input checked="" type="checkbox"/> 8" Castors with central locking	4 pcs
<input checked="" type="checkbox"/> Directional wheel	1 pc
<input checked="" type="checkbox"/> Iv pole previsions	2 pcs
<input checked="" type="checkbox"/> Iv pole	1 set
<input checked="" type="checkbox"/> Central locking break pedal	2 sets
<input checked="" type="checkbox"/> Directional pedal	2 pcs
<input checked="" type="checkbox"/> Transfer mattress	1 pc



Metal material

Baosteel from Fortune Global 500
Panasonic Robotic ensure 360° full
smooth wedding.

**Plastic**

Environmental plastic can be
sterilized in 100°C, tensile upto
30MPa.

**Steels | Painting**


11 Process epoxy painting, ASTM
testing anti-bacterial, paint thickness
0.12 mm, brightness 60°, paint can
resist 50 kg impact.

**Brake**

Break spring keep qualified after
10,000times uninteruped break.

**Hydraulic pumps**

Height adjustment, trendelenburg
and revised tendelenburg through
two USA hydraulic pumps.

**Center wheel**


Can be engaged to facilitate direction
control, enabling users to navigate
the stretcher easily in hectic and
dynamic environments.

**Rising**

Hydraulic treadle for rise.

**Falling**

① Leg section lows down to make
head higher.
② Back section lows down to make
leg section higher.
③ Step on the middle of the pedal,
make the bed base fall down steadily.

**Castors**

Four double side casters of 6 inch
diameter, simultaneously lockable by
the step pedal at both sides.
TPR tire no worn out after running
30KM, save anti-winding hard shell,
united forming without bolts.

**Mattress**


Water proof mattress with 2 belts.

**Side rail**

Bed platform and guardrail are made
of high quality PP material.

**IV pole**

Stainless Steel Column with plastic
hook.





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

Healthward^{uk}

Chengdu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No.1, 5th Floor, Building 5, No.360, Wuke West 5th Road ·
Wuhou District (Chengdu) · China

☎ +86 18021231901 · 🌐 24/7

✉ export@health-ward.com

🌐 www.health-ward.com



Company video



General catalogue